

## Moxi Flow Apoptosis (PE-Annexin V) and Viability (PI) - Staining Protocol

### Reagents:

- Annexin V Binding Buffer ([BioLegend, cat#422201](#))
- Cell Staining Buffer ([BioLegend, cat#420201](#))
- PE-Annexin V conjugate ([BioLegend, cat#640907](#))
- Orflo Flow Reagent ([Orflo Technologies, MXA080](#))
- Moxi Cyte Viability Kit ([Orflo Technologies, MXA055](#))

### Protocol

1. Induce apoptosis in cells by desired method (e.g. 20 $\mu$ M Camptothecin treated, 6+ hours, 37°C) Jurkat cells. Include a control sample of untreated cells.
2. Isolate cells to a single-cell suspension. (If necessary use a protease and/or pipette trituration to break apart the clusters)
3. Wash cells twice in Cell Staining Buffer (2.5mL volume, 300 x g, 5 minutes).
4. Re-suspend pellet to 1 x 10<sup>6</sup> cells/ml in Annexin V Binding Buffer (verify counts with the Moxi Z or Moxi Flow instruments).
5. For each sample, prepare two microcentrifuge tubes. *Note: the tubes can be processed concurrently to reduce overall processing time.*
  - **Tube 1 (PE- Annexin V measurement):**
    - a. Add 100 $\mu$ l of cells
    - b. Add 10  $\mu$ l of PE-Annexin V conjugate. *Notes: Required PE-Annexin V concentration may vary by PE-Annexin manufacturer. We recommend 0.3 – 0.6 $\mu$ g/ml final PE – Annexin V concentration.*
    - c. Gently vortex the cells and incubate for 15 minutes at room temperature (25°C), protected from light.
    - d. Wash cells twice (1mL volume, 300 x g, 5 minutes). with Annexin V Binding Buffer (critical, using a low Ca<sup>2+</sup> or Ca<sup>2+</sup>-free buffer will disrupt Annexin binding)
    - e. Re-suspend pellet in 400 $\mu$ L of Annexin V Binding Buffer
    - f. Add 8  $\mu$ L of Moxi Cyte Flow Reagent. Vortex to mix.
    - g. Run on Moxi Flow using the “Apoptosis (Annexin V-PE)” app within 2 hours of staining, protect from light.
  - **Tube 2 (PI/Viability measurement):**
    - a. Add 60 $\mu$ L of cell suspension to 240 $\mu$ L of Moxi Cyte Viability Kit (5x dilution)
    - b. Gently vortex the cells and incubate for 5 minutes at room temperature (25°C), protected from light.
    - c. Run on Moxi Flow using the “Viability Count (PI)” app within 2 hours of staining, protect from light.
6. As all membrane-permeable (dead) cells will stain positive for PE-Annexin V, the early-stage apoptosis percentage can be calculated as the PE-Annexin V positive percentage (from the first test) minus the PI positive percentage (from the second test).