

Total White Blood Cell (WBC) Count Protocol for use with the Moxi Flow or Zepi Flow

Overview




This approach is a nucleated cell counting method. The technique works by diluting the whole blood, lysing all the cells (including the WBC's), and counting the remaining nuclei. It is designed to provide a rapid means for counting the total WBC counts in a sample. *Note: There is an inherent error in the counts due to contributions of nucleated red blood cells (RBC'S). Ordinarily, this contribution is not significant. However pathological conditions could contribute to higher-than normal nucleated RBC's in blood.*

Reagents

- Zap-oglobin II Lytic Reagent ((Beckman Coulter, [Part # 7546138](#)))
- Ca²⁺-free and Mg²⁺-free PBS or the equivalent
- Anti-coagulant (e.g. Heparin, EDTA, sodium citrate) coated collection tube
- MFS Cassettes ([Orflo Prod #MXC020](#))


Protocol

- Collect >40µL of whole blood into an anti-coagulant coated capillary tube.
- Add exactly 40µL of whole blood to 5mL of PBS/Media.
- Add 3 drops of Zap-oglobin II Lytic Reagent
- Mix the sample several times (10x) through inversion-mixing and allow the sample to sit for a minimum of 3 minutes.
- Turn on the Moxi Flow or Zepi Flow
- Select the appropriate test for the instrument:
 - *Moxi Flow* – Touch “Size Histogram” (red arrow in image to the right)
 - *Zepi Flow* – Touch “Open Flow Cytometry”, select “Histogram” for the Y-Axis column, and touch “Next” (see images below)



Total White Blood Cell (WBC) Count Protocol for use with the Moxi Flow or Zepi Flow

- Load an MFS cassette and close the door(s).
- Prior to loading the sample, mix the sample through slow inversions to ensure proper mixing.
- When the instrument indicates to load the sample (text in the black bar at the top of the screen), open the door and load 75 μ L (minimum) of sample into the cassette loading well.
- Close the door. The test will automatically start.
- Upon test completion, re-scale the axis to 3 – 10 μ m by touching the “3-18 μ m” (magnifying glass icon), and then the “3-10 μ m” (magnifying glass icon, see images below).



- Gate around the nuclei population (typically 4.5 – 10 μ m size range).
- The total, dilution-adjusted counts can be determined by multiplying the Moxi Z concentration by 128 (5.12 mL/.04 mL). Note: This includes ~80 μ L of Zap-oglobin II reagent.