



Introducing 4 New Kits for Moxi Systems

Designed for convenience and affordability, backed by accuracy and precision

MoxiCyte Viability Kit

Propidium Iodide (PI) staining reagent for use with Moxi V and Moxi GO II instruments to provide rapid assessment of cell viability

MoxiCyte Apoptosis Kit

FITC Annexin V and Propidium Iodide (PI) to monitor physiological changes associated with apoptosis

MoxiCyte AOPI Kit

Acridine Orange (AO) to generically label nucleated cells, and Propidium Iodide (PI) to simultaneously determine cellular viability

MoxiCyte Cell Counting Kit

System check beads and protocols to perform highly accurate, volumetric cell counts

Note: All kits include two packs of cassettes (100 tests) and a replacement cassette insert.

www.orflo.com

Scan to order

